

CHROMATOGRAPHY MADE IN GERMANY

PREPARATIVE BIOCOMPATIBLE COLUMN HARDWARE (PEEK)



PREPARATIVE COLUMN HARDWARE: PEEK

Sykam offers a unique range of biocompatible inert PEEK columns of different dimensions.

Pressure Stable up to 250 bar

The PEEK columns are pressure stable and can be used up to 250 bar pressure. This pressure stability makes these columns a perfect alternative to the conventional glass columns.

Axial Compression

The PEEK columns are also available with an axial compression mechanism.



ORDER INFORMATION

Preparative Column Ø10 PEEK

Catalog No	Description	Notes
S009291	PEEK Column Hardware - 10 x 50 mm	
S009292	PEEK Column Hardware - 10 x 75 mm	
S009290	PEEK Column Hardware - 10 x 100 mm	
S009290	PEEK Column Hardware - 10 x 150 mm	
S009290	PEEK Column Hardware - 10 x 200 mm	
S009290	PEEK Column Hardware - 10 x 250 mm	

Catalog No	Description	Notes
S00	Packing Adapter - 10mm	10 x 60 mm

Preparative Column Ø16 PEEK

Catalog No	Description	Notes
S009291	PEEK Column Hardware - 16x 250 mm	
S009292	PEEK Column Hardware - 16 x 500 mm	

Catalog No	Description	Notes
S009291	Packing Adapter - 16mm	16 x 60 mm

Preparative Column Ø20 PEEK

Catalog No	Description	Notes
S009291	PEEK Column Hardware - 20 x 250 mm	
S009292	PEEK Column Hardware - 20 x 500 mm	

Catalog No	Description	Notes
S009291	PEEK Column Hardware - 20 x 250 mm AC	axial compression 240 - 265 mm
S009292	PEEK Column Hardware - 20 x 500 mm AC	axial compression 490 - 515 mm

Catalog No	Description	Notes
S009291	Packing Adapter - 20mm	16 x 60 mm

Preparative Column Ø30 PEEK

Catalog No	Description	Notes
S009291	PEEK Column Hardware - 30 x 250 mm	
S009292	PEEK Column Hardware - 30 x 500 mm	

Catalog No	Description	Notes
S009291	PEEK Column Hardware - 30 x 250 mm AC	axial compression 240 - 265 mm
S009292	PEEK Column Hardware - 30 x 500 mm AC	axial compression 490 - 515 mm

Catalog No	Description	Notes
S00	Packing Adapter - 30mm	30x 60 mm


Preparative Column Ø50 PEEK

Catalog No	Description	Notes
S009291	PEEK Column Hardware - 50 x 250 mm	
S009292	PEEK Column Hardware - 50 x 500 mm	

Catalog No	Description	Notes
S009291	PEEK Column Hardware - 50 x 250 mm AC	axial compression 240 - 265 mm
S009292	PEEK Column Hardware - 50 x 500 mm AC	axial compression 490 - 515 mm

Catalog No	Description	Notes
S00	Packing Adapter - 50mm	50 x 60 mm
S00	Packing Adapter - 50mm AC	50 x 60 mm

- other dimensions are available on request -



Sykam GmbH

Systeme & Komponenten analytischer Meßtechnik

Gewerbering 15
86922 Eresing
Germany

Tel.: ++49 (8193) 93 82 - 0
FAX: ++49 (8193) 93 82 - 20
EMail: info@sykam.com
Web: <http://www.sykam.com>

