

CHROMATOGRAPHY MADE IN GERMANY

NEW!

ION CHROMATOGRAPHY SYSTEM S 155


- ◆ COMPACT ALL-IN-ONE ION CHROMATOGRAPHY SYSTEM
- ◆ CONDUCTIVITY DETECTOR
- ◆ INTEGRATED COLUMN OVEN
- ◆ INTEGRATED ANION AUTO-SUPPRESSOR
- ◆ OPTIONAL VACUUM DEGASSER
- ◆ OPTIONAL ELUENT SELECTION VALVE
- ◆ OPTIONAL INERT GAS SUPPLY REGULATOR



COMPACT IC SYSTEM S 155

The *Compact Ion Chromatograph S 155* is an all-in-one IC system with integrated anion auto-suppressor and column oven.





ORDER INFORMATION

Compact IC System S 155-A

Catalog No	Description	consisting of:
S008544	Compact Ion Chromatography System	S 155 Ion Chromatography System
	for Anion analysis	- integrated Anion auto-suppressor
	without Data System and Separation Column	- integrated column oven
		- integrated manual injection valve S 6120 w. 50 µl sample loop
		- integrated 1-Channel conductivity detector
		- integrated isocratic eluent pump

Compact IC System S 155-C

Catalog No	Description	consisting of:
S008545	Compact Ion Chromatography System	S 155 Ion Chromatography System
	for Cation analysis	- integrated column oven
	without Data System and Separation Column	- integrated manual injection valve S 6120 w. 50 µl sample loop
		- integrated 1-Channel conductivity detector
		- integrated isocratic eluent pump

Upgrade Options

Catalog No	Description	Notes
S003581	Vacuum Degasser	Online 1-Channel vacuum degasser
S008547	Eluent Selection Valve S 6131	Switching valve for 2 eluents and second 2000 ml eluent bottle
S008546	Inert Gas Supply	Gas regulator for inert gas supply and 2000 ml coated eluent bottle

Sykam GmbH

Systeme & Komponenten analytischer Meßtechnik

Gewerbering 15

86922 Eresing

Germany

Tel.: ++49 (8193) 93 82 - 0

FAX: ++49 (8193) 93 82 - 20

E-Mail: info@sykam.com

Web: <http://www.sykam.com>

All specifications are subject to change without notice.